




DOUBLE ACTION PLIERS

INSTRUCTIONS FOR USE

Please read these instructions carefully before application.

Double Action Pliers

GENERAL DESCRIPTION

Double Action Pliers is used for easy and safe removal of FlexiOH-Light Cured Polymer Immobilizer.

INTENDED USE / INDICATION

To be used solely for the removal process of FlexiOH light cured polymer immobilizers.

CONTRAINDICATIONS

None known.

PATIENT POPULATION

Patients undergoing fracture immobilization treatment using FlexiOH - Light Cured Polymer Immobilizer.


INTENDED USERS

A healthcare professional (HCP) such as an orthopedic surgeon, emergency room doctor, physician assistant, orthopedic technician or nurse practitioner.

APPLICATION INSTRUCTION

Follow the steps to remove the FlexiOH Light Cured Polymer Immobilizer before /after the treatment:


1. Cut the fabric liner of the zip using scissors.
2. Select the breaking points as per image.
3. Break the rib of the FlexiOH immobilizer on the opposite side of the zipper line using double action pliers as shown in the image.
4. Apply gentle force to open the immobilizer from the zip side to remove it.



⚡ BREAKING POINT



Keep pliers away from skin



Correct Pliers position

WARNINGS AND PRECAUTIONS

- Keep away from children / Do not allow children to play with Double Action Pliers.
- In case, the box is damaged but the product is intact in the box, the product can be used for the removal of immobilizer.
- Use goggles while Double Action Pliers removing the immobilizer to avoid eye injury.
- Double Action Pliers is for Professional use only.
- To be sold and used only for removal of FlexiOH immobilizer.

ADVERSE EFFECTS

None Known.

STORAGE CONDITIONS

Keep the device in a cool and dry place.

REPORTING OF ADVERSE EVENTS

In case of any event please notify to 'care@orthoheal.com' or to the respective country regulatory authority.

CLEANING GUIDELINE

It is recommended to clean the Double Action Pliers by 70% IPA wipes before every use.

SYMBOL GLOSSARY



Medical Device



Do not use if package is damaged and consult instructions for use



Non-sterile



Consult instructions for use.



Manufacturer



Date of Manufacture



Catalog Number / Reference code



Batch code / Lot Number



Unique device identifier



Authorised representative in the European Union

Copy of the design or technology may infringe IP rights.



FDA
LISTED



MANUFACTURER

JC OrthoHeal Private Limited
P7 - P8/A, Biotech Park Phase - I,
GIDC Savli, Manjusar,
Vadodara - 391 775. Gujarat, India.

Call: +91 85305 30885
Email: india@orthoheal.com

FOR COMPLAINT AND FEEDBACK

Call: +91 911 811 311 0
Email: care@orthoheal.com



AF Pharma Service Europe SL
Muntaner 281,
Barcelona, 08021,
España



2025©Design version 2.0
Revised on 08 January 2025