




FlexiOH[®]

Light Cured Polymer Immobilizer

Absolute washable orthopedic immobilizer for bone fractures and sports injuries.

MEASUREMENT GUIDE (Simplified) in centimeters



SIZE	SHORT ARM	
	Mid Palm Circum.	Forearm Length
3XS	9-11	13-17
2XS		14-18
SXS	11-14	17-20
XS		18-21
SS	14-17	20-23
S		23-26
SM	17-18.5	23-27
M		26-29
SL	18.5-22.5	27-29
L		29-30
SXL	22.5-24.5	29-30
XL		30-33

SIZE	SHORT LEG		FULL ARM		RADIAL GUTTER	ULNAR GUTTER	LONG THUMB SPICA
	Calf Circum.	Head of Metatarsal	Forearm Circum.	Max. Arm Circum.	Wrist Circum.	Wrist Circum.	Wrist Circum.
2XS	18-22	11-13	11.5-15	12-14.5	-	-	-
XS	22-25	13-15	15-18.5	14.5-18	-	-	-
S	25-30	15-18	18.5-23	22-26.5	14-16	14-16	14-16
M	30-40.5	18-22	23-26.5	26.5-31	16-18	16-18	16-18
L	40.5-45	22-25.5	26.5-30	31-36	18-20	18-20	18-20
XL	45-49	25.5-28	-	-	-	-	-

WWW.ORTHOHEAL.COM



Design2026©OrthoHealSCP06

